



DERMAGRIP®

Chemo Protection Gloves

Hand Specific
Straight Fingers
Bisque Finish
Powder Free

Key Features

Beaded cuff

Extra thickness

Chemotherapy tested

Sterilized / Non-sterilized

Undetectable chemical residues

Product Information

Reorder No.	(Sterile) D31XX-19 (Non-Sterile) D31XX-18
Material	Natural rubber latex
Sizes	S-XL
Length (mm)	Min. 290
Target Weight (Size M)	23.3g
Finger Thickness (mm)	Min. 0.46
Colour	Blue
Packing Mode	(Sterile) 100 pr/cs (Non-Sterile) 50 pr/bx, 6 bx/cs

PRODUCT CONFORMANCE

- ASTM D3578
- Meets ASTM D6978 for Chemotherapy testing
- Medical Device: in compliance with 93/42/EEC (CE Class I Sterile)
- EN 455: Part 1/2/3/4
- EN 420:2003+A1:2009, EN ISO 374-1:2016+A1:2018 Type B, EN ISO 374-2:2019, EN ISO 374-4:2019, EN ISO 374-5:2016, EN 16523-1:2015+A1:2018, CE2797
- Health Canada Compliance

QUALITY ASSURANCE

- MDSAP-ISO 13485:2016 Quality Management System
- ISO 9001:2015 Quality Management System
- EN ISO 13485:2016



WRP Asia Pacific Sdn Bhd 1 4 7 8 1 7 V
 Lot 1, Jalan 3, Kawasan Perusahaan Bandar Baru Salak Tinggi,
 43900 Sepang, Selangor Darul Ehsan, MALAYSIA
 Office +60-3-8706-1486
 Email customer_wrp@wrpworld.com

EC REP & Importer in the EU
REMESCO Handelsges.m.b.H.
 Muthgasse 36/19, A-1190 Vienna, AUSTRIA
 Office +43-1-328-50-88
 Fax +43-1-328-29-49
 Email office@remesco.com