



DERMAGRIP®

Examination Gloves *Classic*

NR Latex
Ambidextrous
Powder Free
Fully Textured Surface

Key Features

- Beaded cuff**
- No unpleasant odor**
- Low protein content**
- Strong and durable**
- Excellent fit and feel**
- Double chlorination**

Product Information


Reorder No.	D15xx-10
Material	Natural Rubber Latex
Sizes	XS - XL
Length (mm)	Min. 240
Target Weight (Size M)	7.7g
Finger Thickness (mm)	0.16±0.03
Colour	Natural
Packing Mode	100 pc/bx; 10 bx/cs

PRODUCT CONFORMANCE

- ASTM D3578
- Medical Device: in compliance with MDR (EU)2017/745 (CE Class I)
- EN 455: Part 1/2/3/4
- EN 420:2003+A1:2009, EN ISO 374-1:2016+A1:2018 Type B, EN ISO 374-2:2019, EN ISO 374-4:2019, EN ISO 374-5:2016, EN 16523-1:2015+A1:2018, CE2797
- Health Canada Compliance

QUALITY ASSURANCE

- MDSAP-ISO 13485:2016 Quality Management System
- ISO 9001:2015 Quality Management System
- EN ISO 13485:2016



WRP Asia Pacific Sdn Bhd 1 4 7 8 1 7 V
 Lot 1, Jalan 3, Kawasan Perusahaan Bandar Baru Salak Tinggi,
 43900 Sepang, Selangor Darul Ehsan, MALAYSIA
 Office +60-3-8706-1486
 Email customer_wrp@wrpworld.com

EC REP & Importer in the EU
REMESCO Handelsges.m.b.H.
 Muthgasse 36/19, A-1190 Vienna, AUSTRIA
 Office +43-1-328-50-88
 Fax +43-1-328-29-49
 Email office@remesco.com