



DERMAGRIP®

Nitrile Examination Gloves *Ultra LT*

Nitrile
Ambidextrous
Powder Free
Finger Tip Textured Surface

Key Features

- Beaded cuff**
- Chemotherapy tested**
- No unpleasant odor**
- Puncture resistant**
- Soft and comfortable**
- Superior barrier protection**

Product Information

| | |
|------------------------|--|
| Reorder No. | D15xx-20 & D11xx-23 |
| Material | Nitrile |
| Sizes | XS - XL |
| Length (mm) | Min. 240 |
| Target Weight (Size M) | 3.7g |
| Finger Thickness (mm) | Min. 0.07 |
| Colour | Blue |
| Packing Mode | (D11xx-23) 100pc/bx, 10 bx/cs (D15xx-20) 200pc/bx, 10 bx/cs |

* Size XL is 180 pieces/box

PRODUCT CONFORMANCE

- ASTM D6319
- Meets ASTM D6978 for Chemotherapy testing
- Medical Device: in compliance with MDR (EU)2017/745 (CE Class I)
- EN 455: Part 1/2/3/4
- EN 420:2003+A1:2009, EN ISO 374-1:2016+A1:2018 Type B, EN ISO 374-2:2019, EN ISO 374-4:2019, EN ISO 374-5:2016, EN 16523-1:2015+A1:2018, CE2797
- Food Contact Compliance FDA 21 CFR 177.2600 & EU No. 10/2011
- Health Canada Compliance (D15xx-20)

QUALITY ASSURANCE

- MDSAP-ISO 13485:2016 Quality Management System
- ISO 9001:2015 Quality Management System
- EN ISO 13485:2016



WRP Asia Pacific Sdn Bhd 1 4 7 8 1 7 V
 Lot 1, Jalan 3, Kawasan Perusahaan Bandar Baru Salak Tinggi,
 43900 Sepang, Selangor Darul Ehsan, MALAYSIA
 Office +60-3-8706-1486
 Email customer_wrp@wrpworld.com

EC REP & Importer in the EU
REMESCO Handelsges.m.b.H.
 Muthgasse 36/19, A-1190 Vienna, AUSTRIA
 Office +43-1-328-50-88
 Fax +43-1-328-29-49
 Email office@remesco.com