



# DERMAGRIP®

NITRILE EXAMINATION GLOVES **ULTRA PLUS**

**Nitrile**  
**Ambidextrous**  
**Powder Free**  
**Finger Textured Surface**

## Key Features

- Beaded cuff**
- Puncture resistant**
- Low chemical residue**
- Chemotherapy tested**
- Enabling good dexterity**
- Longer cuff for extra protection**

## Product Information

|                        |  |
|------------------------|--|
| Reorder No.            | D14xx-27 & D14xx-28  |
| Material               | Nitrile  |
| Sizes                  | XS - XL  |
| Length (mm)            | Min. 270   |
| Target Weight (Size M) | 4.4g   |
| Finger Thickness (mm)  | 0.11 ± 0.02  |
| Colour                 | Blue   |
| Packing Mode           | (D14xx-27) 100 pc/bx, 10 bx/cs<br>(D14xx-28) 200 pc/bx, 10 bx/cs |

\* Size XL is 180 pieces/box

## PRODUCT CONFORMANCE

- ASTM D6319
- Meets ASTM D6978 for Chemotherapy testing
- Medical Device: in compliance with MDR (EU)2017/745 (CE Class I)
- EN 455: Part 1/2/3/4
- EN 420:2003+A1:2009, EN ISO 374-1:2016+A1:2018 Type B, EN ISO 374-2:2019, EN ISO 374-4:2019, EN ISO 374-5:2016, EN 16523-1:2015+A1:2018, CE2797
- Food Contact Compliance FDA 21 CFR 177.2600 & EU No. 10/2011

## QUALITY ASSURANCE

- MDSAP-ISO 13485:2016 Quality Management System
- ISO 9001:2015 Quality Management System
- EN ISO 13485:2016



**WRP Asia Pacific Sdn Bhd** 1 4 7 8 1 7 V  
 Lot 1, Jalan 3, Kawasan Perusahaan Bandar Baru Salak Tinggi,  
 43900 Sepang, Selangor Darul Ehsan, MALAYSIA  
 Office +60-3-8706-1486  
 Email customer\_wrp@wrpworld.com

**EC REP & Importer in the EU**  
**REMESCO Handelsges.m.b.H.**  
 Muthgasse 36/19, A-1190 Vienna, AUSTRIA  
 Office +43-1-328-50-88  
 Fax +43-1-328-29-49  
 Email office@remesco.com