



** 3D Diamond Texture

E-DUO™

Multipurpose Gloves Extended Cuff

Nitrile

Ambidextrous

Powder Free

Diamond Textured Surface

Key Features

Beaded cuff

Chemical resistant

Strong and durable

Chemotherapy tested

Embossed diamond texture

Excellent dexterity and elasticity

Product Information

Reorder No.	E105xx-xx
Material	Nitrile
Sizes	S - XXL
Length (mm)	Min. 290
Target Weight (Size M)	10.5g
Finger Thickness (mm)	0.27 + 0.03
Colour	Blue, Green, Black, Orange
Packing Mode	50 pc/bx, 10 bx/cs

PRODUCT CONFORMANCE

- ASTM D6319
- Meets ASTM D6978 for Chemotherapy testing
- Medical Device: in compliance with MDR (EU)2017/745 (CE Class I)
- EN 455: Part 1/2/3/4
- EN 420:2003+A1:2009, EN ISO 374-1:2016+A1:2018 Type B, EN ISO 374-2:2019, EN ISO 374-4:2019, EN ISO 374-5:2016, EN 16523-1:2015+A1:2018, CE2797
- Food Contact Compliance FDA 21 CFR 177.2600 & EU No. 10/2011

QUALITY ASSURANCE

- MDSAP-ISO 13485:2016 Quality Management System
- ISO 9001:2015 Quality Management System
- EN ISO 13485:2016



WRP Asia Pacific Sdn Bhd 1 4 7 8 1 7 V
 Lot 1, Jalan 3, Kawasan Perusahaan Bandar Baru Salak Tinggi,
 43900 Sepang, Selangor Darul Ehsan, MALAYSIA
 Office +60-3-8706-1486
 Email customer_wrp@wrpworld.com

EC REP & Importer in the EU
REMESCO Handelsges.m.b.H.
 Muthgasse 36/19, A-1190 Vienna, AUSTRIA
 Office +43-1-328-50-88
 Fax +43-1-328-29-49
 Email office@remesco.com