



E-DUO™

Multipurpose Gloves **PLUS**

Nitrile
Ambidextrous
Powder Free
Diamond Textured Surface

Key Features

Beaded cuff

Chemical resistant

Strong and durable

Chemotherapy tested

Embossed diamond texture

Excellent dexterity and elasticity

Product Information

Reorder No.	E85xx-xx
Material	Nitrile
Sizes	S - XXL
Length (mm)	Min. 240
Target Weight (Size M)	8.5g
Finger Thickness (mm)	Min 0.23 ± 0.03
Colour	Blue, Green, Black, Orange
Packing Mode	50 pc/bx, 10 bx/cs

PRODUCT CONFORMANCE

- ASTM D6319
- Meets ASTM D6978 for Chemotherapy testing
- Medical Device: in compliance with MDR (EU)2017/745 (CE Class I)
- EN 455: Part 1/2/3/4
- EN 420:2003+A1:2009, EN ISO 374-1:2016+A1:2018 Type B, EN ISO 374-2:2019, EN ISO 374-4:2019, EN ISO 374-5:2016, EN 16523-1:2015+A1:2018, CE2797
- Food Contact Compliance FDA 21 CFR 177.2600 & EU No. 10/2011

QUALITY ASSURANCE

- MDSAP-ISO 13485:2016 Quality Management System
- ISO 9001:2015 Quality Management System
- EN ISO 13485:2016



WRP Asia Pacific Sdn Bhd 1 4 7 8 1 7 V
 Lot 1, Jalan 3, Kawasan Perusahaan Bandar Baru Salak Tinggi,
 43900 Sepang, Selangor Darul Ehsan, MALAYSIA
 Office +60-3-8706-1486
 Email customer_wrp@wrpworld.com

EC REP & Importer in the EU
REMESCO Handelsges.m.b.H.
 Muthgasse 36/19, A-1190 Vienna, AUSTRIA
 Office +43-1-328-50-88
 Fax +43-1-328-29-49
 Email office@remesco.com