



PROFEEL®

DOUBLE GLOVING SYSTEM

FOR PUNCTURE DETECTION

DHD™ POLYISOPRENE SURGICAL UNDERGLOVE AND OVERGLOVE
NON - LATEX POWDER FREE

Hand Specific
Powder Free
Straight Fingers
Damp-Hand-Donnable (DHD™)
Micro Roughened Surface

Key Features

Sterilized

Latex free

Smooth finish

Polymer coated

Chemotherapy tested

Standard thickness / Micro thickness

Product Information

	<i>Overglove</i>	<i>Underglove</i>
Reorder No.	P33xx-62	P33xx-62
Material	Polyisoprene	Polyisoprene
Sizes	5.5 - 9.0	5.5 - 9.0
Length (mm)	± 290	± 290
Target Weight (Size 7.5)	25.91g	25.91g
Finger Thickness (mm)	0.24 ± 0.02	0.20 ± 0.02
Colour	Natural	Green
Packing Mode	25pr/bx; 4 bx/cs	25pr/bx; 4 bx/cs

PRODUCT CONFORMANCE

- ASTM D3577
- Meets ASTM D6978 for Chemotherapy testing
- Medical Device: in compliance with 93/42/EEC
- EN 455: Part 1/2/3/4
- EN 420:2003+A1:2009, EN ISO 374-1:2016+A1:2018 Type B, EN ISO 374-2:2019, EN ISO 374-4:2019, EN ISO 374-5:2016, EN 16523-1:2015+A1:2018, CE2797

QUALITY ASSURANCE

- MDSAP-ISO 13485:2016 Quality Management System
- ISO 9001:2015 Quality Management System
- EN ISO 13485:2016



WRP Asia Pacific Sdn Bhd 1 4 7 8 1 7 V
Lot 1, Jalan 3, Kawasan Perusahaan Bandar Baru Salak Tinggi,
43900 Sepang, Selangor Darul Ehsan, MALAYSIA
Office +60-3-8706-1486
Email customer_wrp@wrpworld.com

EC REP & Importer in the EU
REMESCO Handelsges.m.b.H.
Muthgasse 36/19, A-1190 Vienna, AUSTRIA
Office +43-1-328-50-88
Fax +43-1-328-29-49
Email office@remesco.com