



**DERMAGRIP**



**NR LATEX EXAMINATION GLOVES,**

**POWDER FREE, NON-STERILE**

Reorder #	Size: XS S M L XL	D1500-50 D1501-50 D1502-50 D1503-50 D1504-50												
Product Part #	156xx.10021746													
510K # / MDL #	N/A													
Design and Feature	Ambidextrous Fingers-Textured E5,240mm, Strip and Pack													
Type	Powder free and non-sterile examination gloves													
Material	Natural rubber latex													
Color	Natural (i.e. off white to yellow)													
Surface Treatment	Strip & Pack													
Protein Content	This latex glove contains 50 microgram or less of total extractable protein per gram													
Powder Free Residue (mg/glove)	≤ 2													
Dimension	Size: XS S M L XL	<table border="0"> <tr> <td><u>Palm Width (mm)</u></td> <td><u>Length (mm)</u></td> </tr> <tr> <td>≤80</td> <td>240 min</td> </tr> <tr> <td>80 ± 10</td> <td>240 min</td> </tr> <tr> <td>95 ± 10</td> <td>240 min</td> </tr> <tr> <td>110 ± 10</td> <td>240 min</td> </tr> <tr> <td>≥110</td> <td>240 min</td> </tr> </table>	<u>Palm Width (mm)</u>	<u>Length (mm)</u>	≤80	240 min	80 ± 10	240 min	95 ± 10	240 min	110 ± 10	240 min	≥110	240 min
<u>Palm Width (mm)</u>	<u>Length (mm)</u>													
≤80	240 min													
80 ± 10	240 min													
95 ± 10	240 min													
110 ± 10	240 min													
≥110	240 min													
Single-wall Thickness (mm) *All sizes	Finger Palm Cuff	0.14 ± 0.03 0.12 ± 0.03 0.10 ± 0.02												
"Typical" Physical Properties	Tensile Strength (MPa) Ultimate Elongation (%) Stress at 500% (MPa) Force at Break (N)	<table border="0"> <tr> <td><u>Before Aging</u></td> <td><u>After Aging</u></td> </tr> <tr> <td>Min. 18</td> <td>Min. 14</td> </tr> <tr> <td>Min. 650</td> <td>Min. 500</td> </tr> <tr> <td>Max. 5.5</td> <td>Not applicable</td> </tr> <tr> <td>≥ 6</td> <td>≥ 6</td> </tr> </table>	<u>Before Aging</u>	<u>After Aging</u>	Min. 18	Min. 14	Min. 650	Min. 500	Max. 5.5	Not applicable	≥ 6	≥ 6		
<u>Before Aging</u>	<u>After Aging</u>													
Min. 18	Min. 14													
Min. 650	Min. 500													
Max. 5.5	Not applicable													
≥ 6	≥ 6													
Packing Mode	100 gloves (by weight per dispenser) 10 dispensers per carton /shipping case													
Glove Marking	No marking.													
Lot # identification on Finished Goods	Lot Number Structure: YMMPPPPSS (9 digits) Y = Year of Packing MM = Month of Packing PPPP = WRP's Packing Work Order # (PWO) SS = Size (00=XS, 01=S, 02=M, 03=L and 04=XL)													
Recommended Shelf Life	3 years upon manufactured date.													
Pre-shipment Inspection  *Single-Normal Sampling Plan	Dimension Physical Properties 1000ml Water Leak Protein Content Powder Free Residue Major Visual Inspection Minor Visual Inspection	N=13, Median; S-2 AQL 4.0 N=13, Median; S-2 AQL 4.0 G-I, AQL 1.5 N=8 N=5 G-I, AQL 2.5 G-I, AQL 4.0												
Product Conformance	<ul style="list-style-type: none"> <li>Personal Protective Equipment of Complex Design Category III, in compliance with European Regulation (EU) 2016/425, tested to EN 420:2003+A1:2009, EN ISO 374-1:2016 Type B, EN 374-2:2014, EN 16523-1:2015, EN 374-4:2013 &amp; EN ISO 374-5:2016, CE2797; CE Class I Medical Device, in compliance with MDR (EU) 2017/745, type tested to EN 455-1:2000, EN 455-2:2015, EN 455-3:2015, EN 455-4:2009</li> <li>ASTM D3578</li> </ul>													
Quality Assurance	<ul style="list-style-type: none"> <li>US FDA Quality System Regulation (QSR)</li> <li>ISO9001:2015 Quality Management Systems</li> <li>EN ISO13485:2016 &amp; ISO13485:2016 Quality Management Systems</li> </ul>													

Note: The above information is a guideline of typical performance values and characteristic of the product; and not to be used as actual product specifications.