



SPECIFICATION SHEET

DERMAGRIP®

Nitrile Examination Gloves, Powder Free, Non-Sterile

Size	Reorder #	Dimension		
		Target weight (g/pcs)	Palm Width (mm)	Length (mm)
Extra Small	D1500-17	4.1 ± 0.3	≤ 80	Min. 240
Small	D1501-17	4.8 ± 0.3	80 ± 10	Min. 240
Medium	D1502-17	5.5 ± 0.3	95 ± 10	Min. 240
Large	D1503-17	6.2 ± 0.3	110 ± 10	Min. 240
Extra Large	D1504-17	6.9 ± 0.3	≥ 110	Min. 240
		<i>*non-release criteria</i>		
Product Part #		142xx.10231043		
Design and Feature		Ambidextrous, textured surface at fingers and beaded cuff		
Type		Powder free and non-sterile examination gloves		
Material		100% Nitrile (Acrylonitrile-butadiene)		
Colour		Blue		
Surface Treatment		On-line chlorination on donning side		
Powder Residue (mg/glove)		≤ 2		
Single-wall Thickness (mm) <i>*All sizes</i>	Finger Palm Cuff	Specification Min. 0.14 Min. 0.08 Min. 0.08		
Physical Properties	Tensile Strength (MPa) Ultimate Elongation (%) Force at break (N)	Unaged Min. 14 Min. 500 ≥ 9	Aged Min. 14 Min. 400 ≥ 6	
Packing Mode		100 gloves by weight per dispenser, 10 dispensers per carton		
Glove Marking		No marking		
Lot # Identification on Finished Goods		Lot Number Structure : YMMPPPPSS (9 Digits) Y = Year of packing (e.g. 3 for year 2023, 4 for year 2024, etc.) MM = Month of packing (e.g. 09 for Sept., 10 for Oct., etc.) PPPP = WRP's Packing Work Order# (PWO) SS = Size (00=XS, 01=S, 02=M, 03=L, 04=XL and 05=XXL)		
Shelf Life Claim		5 years upon manufacturing date		
Pre-shipment Inspection <i>*Single-Normal Sampling Plan</i>	Dimension Physical Properties 1000ml Water Leak Powder Residue Major Visual Inspection Minor Visual Inspection	N=13, Median; S-2, AQL 4.0 N=13, Median; S-2, AQL 4.0 G-1, AQL 1.5 N=5 G-1, AQL 2.5 G-1, AQL 4.0		
Product Conformance		<ul style="list-style-type: none"> Medical Device: in compliance with MDR (EU)2017/745 (CE Class I) Personal Protective Equipment of Complex Design Category III, in compliance with European Regulation (EU) 2016/425, Regulation 2016/425 on personal protective equipment, as amended to apply in GB, tested to EN ISO 21420:2020, EN ISO 374-1:2016+A1:2018 Type B, EN ISO 374-2:2019, EN ISO 374-4:2019, EN ISO 374-5:2016, EN 16523-1:2015+A1:2018, CE2797 EN 455 Part 1/2/3/4 ASTM D6319, ASTM D6978 (Chemotherapy testing) Food Safe Compliance (FDA 21 CFR 177.2600 & EU No. 10/2011) Health Canada Compliance 		
Quality Assurance		<ul style="list-style-type: none"> US FDA Quality System Regulation (QSR) MDSAP ISO 13485:2016 Quality Management System ISO 9001:2015 Quality Management System EN ISO 13485:2016 ISO 14001:2015, Environmental Management System 		

Note: The above information is a guideline of typical performance values and characteristics of the product.

WRP ASIA PACIFIC SDN. BHD.

 

Name : Rosnaini Binti Rashidin
Title : Product Management Manager