



SPECIFICATION SHEET

**DERMAGRIP<sup>®</sup> ULTRA PLUS**

**NITRILE EXAMINATION GLOVES, POWDER FREE, NON-STERILE**

Size	Reorder #	Dimension		
		Target weight (g/pcs)	Palm Width (mm)	Length (mm)
Extra Small	D1400-28	3.7 ± 0.3	≤ 80	Min. 270
Small	D1401-28	4.0 ± 0.3	80 ± 10	Min. 270
Medium	D1402-28	4.4 ± 0.3	95 ± 10	Min. 270
Large	D1403-28	5.2 ± 0.3	110 ± 10	Min. 270
Extra Large	D1404-28	5.5 ± 0.3	≥ 110	Min. 270
		<i>*non-release criteria</i>		
Product Part No.		142xx.102511030		
Design and Feature		Ambidextrous, finger textured and beaded cuff		
Type		Powder free and non-sterile examination gloves		
Material		100% nitrile (Acrylonitrile-butadiene)		
Colour		Blue		
Surface Treatment		Chlorination on donning side		
Powder Residue (mg/glove)		≤ 2		
Single-wall Thickness (mm)	Finger	<u>Specification</u>		
*All sizes	Palm	0.11 ± 0.02		
	Cuff	0.08 ± 0.02		
		0.07 ± 0.02		
Physical Properties	Tensile Strength (MPa)	Unaged	Aged	
		Min. 14	Min. 14	
Ultimate Elongation (%)	Force at Break (N)	Min. 500	Min. 400	
		≥ 6	≥ 6	
Packing Mode		200 gloves by weight per dispenser and 10 dispensers per carton (XS-L) 180 gloves by weight per dispenser and 10 dispensers per carton (XL)		
Glove Marking		No marking		
Lot # Identification on Finished Goods		Lot Number Structure : YMMPPPPSS (9 Digits) Y = Year of packing (e.g. 3 for year 2023, 4 for year 2024, etc.) MM = Month of packing (e.g. 09 for Sept., 10 for Oct., etc.) PPPP = WRP's Packing Work Order# (PWO) SS = Size (00=XS, 01=S, 02=M, 03=L and 04=XL)		
Shelf Life Claim		5 years upon manufacturing date		
Pre-shipment Inspection	Dimension	N=13, Median; S-2, AQL 4.0		
	Physical Properties	N=13, Median; S-2, AQL 4.0		
*Single-Normal Sampling Plan	1000ml Water Leak	G-I, AQL 1.5		
	Protein Content	N=8 (EN 455); N=3 (ASTM D3578)		
	Powder Residue	N=5		
	Major Visual Inspection	G-I, AQL 2.5		
Minor Visual Inspection	G-I, AQL 4.0			
Product Conformance		<ul style="list-style-type: none"> <li>Medical Device: in compliance with MDR (EU)2017/745 (CE Class I)</li> <li>EN 420:2003+A1:2009, EN ISO 374-1:2016+A1:2018 Type B, EN ISO 374-2:2019, EN ISO 374-4:2019, EN ISO 374-5:2016, EN 16523-1:2015+A1:2018, CE2797</li> <li>EN 455 Part 1/2/3/4</li> <li>ASTM D6319</li> <li>Food Safe Compliance (FDA 21 CFR 177.2600 &amp; EU No. 10/2011)</li> </ul>		
Quality Assurance		<ul style="list-style-type: none"> <li>MDSAP ISO 13485:2016 Quality Management System</li> <li>ISO 9001:2015 Quality Management System</li> <li>EN ISO 13485:2016</li> <li>ISO 14001:2015, Environmental Management System</li> </ul>		

Note: The above information is a guideline of typical performance values and characteristic of the product.

WRP ASIA PACIFIC SDN. BHD.

Name : Rosnaini Binti Rashidin  
Title : Product Management Manager