




**DERMAGRIP®**  
HIGH RISK EXAMINATION GLOVES  
NR LATEX POWDER FREE

**NR Latex**  
**Ambidextrous**  
**Powder Free**  
**Finger Textured Surface**

**Key Features**

- Beaded cuff**
- No unpleasant odor**
- Chemotherapy tested**
- Extra-long and non-sterile protection**
- Thicker and stronger**
- Protection in high risk situations**

**Product Information**

Reorder No.	D14xx-14
Material	Natural Rubber Latex
Sizes	S - XL
Length (mm)	± 290
Target Weight (Size M)	16g
Finger Thickness (mm)	0.36 ± 0.04
Colour	Blue
Packing Mode	50 pc/bx, 10 bx/cs

**PRODUCT CONFORMANCE**

- ASTM D3578
- Meets ASTM D6978 for Chemotherapy testing
- Medical Device: in compliance with MDR (EU)2017/745 (CE Class I)
- EN 455: Part 1/2/3/4
- EN ISO 21420:2020, EN ISO 374-1:2016+A1:2018 Type A, EN ISO 374-2:2019, EN ISO 374-4:2019, EN ISO 374-5:2016, EN 16523-1:2015+A1:2018, ISO 16604:2004, CE2797
- Food Contact Compliance EU No. 10/2011
- Health Canada Compliance

**QUALITY ASSURANCE**

- MDSAP-ISO 13485:2016 Quality Management System
- ISO 9001:2015 Quality Management System
- EN ISO 13485:2016
- ISO 14001:2015 Environmental Management System



**WRP Asia Pacific Sdn Bhd**

Lot 1, Jalan 3, Kawasan Perusahaan Bandar Baru Salak Tinggi  
43900 Sepang, Selangor Darul Ehsan, Malaysia  
Phone +60-3-8706-1486  
Email [inquiry@wrpworld.com](mailto:inquiry@wrpworld.com)

**EC REP** & Importer in the EU

Remesco Handelsges.m.b.H.  
Muthgasse 36/19, A-1190 Vienna, Austria  
Phone +43-1-328-50-88  
Email [office@remesco.com](mailto:office@remesco.com)