



# DERMAGRIP®

Nitrile Examination Gloves **Ultra LS**

**Nitrile**  
**Ambidextrous**  
**Powder Free**  
**Finger Tip Textured Surface**

## Key Features

- Beaded cuff**
- Puncture resistant**
- No unpleasant odor**
- Chemotherapy tested**
- Soft and comfortable**

## Product Information

Reorder No.	D11xx-24 & D11xx-25
Material	Nitrile
Sizes	XS - XL
Length (mm)	Min. 240
Target Weight (Size M)	3.5g
Finger Thickness (mm)	0.08 ± 0.03
Colour	Blue
Packing Mode	(D11xx-25) 100 pc/bx, 10 bx/cs (D11xx-24) 200 pc/bx, 10 bx/cs

\* Size XL is 180 pieces/box

## PRODUCT CONFORMANCE

- ASTM D6319
- Meets ASTM D6978 for Chemotherapy testing
- Medical Device: in compliance with MDR (EU)2017/745 (CE Class I)
- EN 455: Part 1/2/3/4
- EN ISO 21420:2020, EN ISO 374-1:2016+A1:2018 Type B, EN ISO 374-2:2019, EN ISO 374-4:2019, EN ISO 374-5:2016, EN 16523-1:2015+A1:2018, ISO 16604:2004, CE2797
- Food Contact Compliance FDA 21 CFR 177.2600 & EU No. 10/2011
- Health Canada Compliance

## QUALITY ASSURANCE

- MDSAP-ISO 13485:2016 Quality Management System
- ISO 9001:2015 Quality Management System
- EN ISO 13485:2016
- ISO 14001:2015 Environmental Management System



## WRP Asia Pacific Sdn Bhd

Lot 1, Jalan 3, Kawasan Perusahaan Bandar Baru Salak Tinggi  
 43900 Sepang, Selangor Darul Ehsan, Malaysia  
 Phone +60-3-8706-1486  
 Email [inquiry@wrpworld.com](mailto:inquiry@wrpworld.com)

**EC REP** & Importer in the EU

Remesco Handelsges.m.b.H.  
 Muthgasse 36/19, A-1190 Vienna, Austria  
 Phone +43-1-328-50-88  
 Email [office@remesco.com](mailto:office@remesco.com)