



DERMAGRIP®

Examination Gloves **Ultra**

NR Latex
Ambidextrous
Powder Free
Finger Tip Textured Surface

Key Features

- Beaded cuff**
- Superior strength**
- Puncture resistant**
- Single chlorination**
- Excellent fit and feel**

Product Information

Reorder No.	D10xx-20
Material	Natural Rubber Latex
Sizes	XS - XL
Length (mm)	Min. 240
Target Weight (Size M)	6.5g
Finger Thickness (mm)	0.14 ± 0.03
Colour	Natural
Packing Mode	100 pc/bx; 10 bx/cs

PRODUCT CONFORMANCE

- ASTM D3578
- Medical Device: in compliance with MDR (EU)2017/745 (CE Class I)
- EN 455: Part 1/2/3/4
- EN ISO 21420:2020, EN ISO 374-1:2016+A1:2018 Type B, EN ISO 374-2:2019, EN ISO 374-4:2019, EN ISO 374-5:2016, EN 16523-1:2015+A1:2018, ISO 16604:2004, CE2797
- Health Canada Compliance

QUALITY ASSURANCE

- MDSAP-ISO 13485:2016 Quality Management System
- ISO 9001:2015 Quality Management System
- EN ISO 13485:2016
- ISO 14001:2015 Environmental Management System



WRP Asia Pacific Sdn Bhd

Lot 1, Jalan 3, Kawasan Perusahaan Bandar Baru Salak Tinggi
 43900 Sepang, Selangor Darul Ehsan, Malaysia
 Phone +60-3-8706-1486
 Email inquiry@wrpworld.com

EC REP & Importer in the EU

Remesco Handelsges.m.b.H.
 Muthgasse 36/19, A-1190 Vienna, Austria
 Phone +43-1-328-50-88
 Email office@remesco.com