



PROFEEL®

**DHD™ Polyisoprene Green Powder Free
Non-Latex Surgical Gloves**

**Damp Hand Donnable (DHD™)
Powder-Free
Excellent grip
Exceptional Comfort**

Key Features

- Latex Free**
- Polymer Coated**
- Micro-roughened surface**
- Optimal instrument handling and grip**
- Beaded Cuff**
- Sterile**

Product Information

Reorder No.	P33xx-56G
Material	Polyisoprene
Sizes	5.5 - 9.0
Length (mm)	± 290
Target Weight (Size 7.5)	11.2g
Finger Thickness (mm)	0.22 ± 0.03
Colour	PMS 325 Green
Packing Mode	50 pr/bx; 4 bx/cs

PRODUCT CONFORMANCE

- ASTM D3577
- Medical Device: in compliance with 93/42/EEC
- EN 455: Part 1/2/3/4
- Health Canada Compliance

QUALITY ASSURANCE

- MDSAP ISO 13485:2016 Quality Management System
- ISO 9001:2015 Quality Management System
- EN ISO 13485:2016
- ISO 14001:2015 Environmental Management System



WRP Asia Pacific Sdn Bhd

Lot 1, Jalan 3, Kawasan Perusahaan Bandar Baru Salak Tinggi
43900 Sepang, Selangor Darul Ehsan, Malaysia
Phone +60-3-8706-1486
Email inquiry@wrpworld.com

EC REP & Importer in the EU

Remesco Handelsges.m.b.H.
Muthgasse 36/19, A-1190 Vienna, Austria
Phone +43-1-328-50-88
Email office@remesco.com