



DERMAGRIP®

NR LATEX EXAMINATION GLOVES, POWDER FREE, STERILE

Reorder #	Size: : XS S M L XL	Pair-Pack D1500-13 D1501-13 D1502-13 D1503-13 D1504-13												
Product Part #	145xx.10022745													
510K # / MDL #	N/A													
Design and Feature	Ambidextrous, finger-textured surface (E4) and beaded cuff													
Type	Powder free and sterile examination gloves													
Material	Natural rubber latex													
Color	Natural (i.e. off white to yellow)													
Surface Treatment	Light chlorination on both donning and application surfaces													
Protein Content	This latex glove contains 50 microgram or less of total extractable protein per gram													
Powder Free Residue (mg/glove)	≤ 2													
Dimension	Size: : XS S M L XL	<table border="1"> <thead> <tr> <th><u>Palm Width</u></th> <th><u>Length</u></th> </tr> </thead> <tbody> <tr> <td>80 ± 10</td> <td>Min. 240</td> </tr> <tr> <td>85 ± 10</td> <td>Min. 240</td> </tr> <tr> <td>95 ± 10</td> <td>Min. 240</td> </tr> <tr> <td>110 ± 10</td> <td>Min. 240</td> </tr> <tr> <td>115 ± 10</td> <td>Min. 240</td> </tr> </tbody> </table>	<u>Palm Width</u>	<u>Length</u>	80 ± 10	Min. 240	85 ± 10	Min. 240	95 ± 10	Min. 240	110 ± 10	Min. 240	115 ± 10	Min. 240
<u>Palm Width</u>	<u>Length</u>													
80 ± 10	Min. 240													
85 ± 10	Min. 240													
95 ± 10	Min. 240													
110 ± 10	Min. 240													
115 ± 10	Min. 240													
Single-wall Thickness (mm) *All sizes	Finger Palm Cuff	<table border="1"> <thead> <tr> <th><u>Spec</u></th> </tr> </thead> <tbody> <tr> <td>0.16 ± 0.03</td> </tr> <tr> <td>0.14 ± 0.03</td> </tr> <tr> <td>0.11 ± 0.03</td> </tr> </tbody> </table>	<u>Spec</u>	0.16 ± 0.03	0.14 ± 0.03	0.11 ± 0.03								
<u>Spec</u>														
0.16 ± 0.03														
0.14 ± 0.03														
0.11 ± 0.03														
Physical Properties	Tensile Strength (MPa) Ultimate Elongation (%) Stress at 500% (MPa) Force at Break (N)	<table border="1"> <thead> <tr> <th><u>Before Aging</u></th> <th><u>After Aging</u></th> </tr> </thead> <tbody> <tr> <td>Min. 18</td> <td>Min. 14</td> </tr> <tr> <td>Min. 650</td> <td>Min. 500</td> </tr> <tr> <td>Max 5.5</td> <td>Not applicable</td> </tr> <tr> <td>≥ 6</td> <td>≥ 6</td> </tr> </tbody> </table>	<u>Before Aging</u>	<u>After Aging</u>	Min. 18	Min. 14	Min. 650	Min. 500	Max 5.5	Not applicable	≥ 6	≥ 6		
<u>Before Aging</u>	<u>After Aging</u>													
Min. 18	Min. 14													
Min. 650	Min. 500													
Max 5.5	Not applicable													
≥ 6	≥ 6													
Packing Mode	2 gloves per inner wrap, 1 inner wrap per pouch, 50 pouches per dispenser and 6 dispensers per carton													
Glove Marking	No marking.													
Lot # identification on Finished Goods	Lot Number Structure: YMMPPPPSS (9 digits) Y = Year of Packing MM = Month of Packing PPPP = WRP's Packing Work Order # (PWO) SS = Size (00=XS, 01=S, 02=M, 03=L and 04=XL)													
Recommended Shelf Life	3 years upon manufactured date.													
Pre-shipment Inspection *Single-Normal Sampling Plan	Dimension Physical Properties 1000ml Water Leak Protein Content Powder Free Residue Major Visual Inspection Minor Visual Inspection	N=13, Median N=13, Median G-I, AQL 1.5 N=3 N=5 G-I, AQL 2.5 G-I, AQL 4.0												
Product Conformance	<ul style="list-style-type: none"> Medical Device: In compliance with European MDR (EU) 2017/745 (CE CLASS Is) EN455 Parts 1, 2, 3 and 4 Personal Protective Equipment of Complex Design Category III, in European Regulation (EU) 2016/425, tested to EN 420:2003+A1:2009, EN ISO 374-1:2016 Type B, EN 374-2:2014, EN 374-4:2013, EN ISO 374-5:2016 & EN 16523-1:2015, CE 2797 ASTM D3578 													
Quality Assurance	<ul style="list-style-type: none"> US FDA Quality System Regulation (QSR) ISO9001:2015 Quality Management Systems EN ISO13485:2016 & ISO13485:2016 Quality Management Systems 													

Note: The above information is a guideline of typical performance values and characteristic of the product; and not to be used as actual product specifications.