



WRP Asia Pacific Sdn Bhd

14 7 8 17 V

Lot 1, Jalan 3, Kawasan Perusahaan

Bandar Baru Salak Tinggi,

43900 Sepang,

Selangor Darul Ehsan, MALAYSIA

Office +60-3-8706 1486

Facsimile +60-3-8706 1557

Email customer_wrp@wrpworld.com

Website www.wrpworld.com

DECLARATION OF FOOD CONTACT COMPLIANCE

Regulation : FDA EXTRACTION TEST: CFR 177.2600
: Code of Federal Regulation, Title 21
: Chapter 1 – Food and Drug Administration
: Part 177.2600 – Rubber articles intended for repeated use

Type of Glove : Black Nitrile Powder Free Examination Gloves

Validity : Five (5) Years

Date of Analysis : 22nd September 2020

OVERALL MIGRATION

METHOD OF TEST:

The sample was analysed for the following tests according to FDA 21 CFR 177.2600 (e) & (f).

1. Preparation of Test Specimen

The sample was cut and tested as received.

2. Total Extractives with Food Stimulants (Ultrapure Water, n-Hexane):

The sample was extracted with the food stimulants and extracted at the reflux temperature for 7 hours and then for another further succeeding 2 hours. The extract was then evaporated to dryness and the total extractives were then measures using an analytical balance.



WRP Asia Pacific Sdn Bhd

147817V

Lot 1, Jalan 3, Kawasan Perusahaan
Bandar Baru Salak Tinggi,
43900 Sepang,
Selangor Darul Ehsan, MALAYSIA

Office +60-3-8706 1486

Facsimile +60-3-8706 1557

Email customer_wrp@wrpworld.com

Website www.wrpworld.com

RESULTS:

EXTRACTION TIMES	SOLVENT	REQUIREMENT (mg/in.²) max.	INFERRED RESULT
First 7 hours	Distilled Water	20	PASS
Next 2 hours	Distilled Water	1	PASS
First 7 hours	n-Hexane	175	PASS
Next 2 hours	n-Hexane	4	PASS

End of report

Reference Report : Test Report No.: PN 155881

Test Facility : Akron Rubber Development Laboratory, Inc (ARDL)

Prepared by,

Approved by,


.....
Rosnaini Rashidin
Product Management Manager


.....
Angie Ng Geok Chin
Research & Development Manager