




**SPECIFICATION SHEET**



**PCP NR LATEX POWDER FREE SURGICAL GLOVES, STERILE**

Size	Reorder#	Dimension		
		Target weight (g/pcs)	Palm Width (mm)	Length (mm)
5.5	P3155-21	9.2 ± 0.5	72 ± 4	min. 278
6.0	P3160-21	9.8 ± 0.5	77 ± 5	min. 280
6.5	P3165-21	10.3 ± 0.5	83 ± 5	min. 280
7.0	P3170-21	10.9 ± 0.5	89 ± 5	min. 283
7.5	P3175-21	11.5 ± 0.5	95 ± 5	min. 287
8.0	P3180-21	12.1 ± 0.5	102 ± 6	min. 288
8.5	P3185-21	12.7 ± 0.5	108 ± 6	min. 290
9.0	P3190-21	13.2 ± 0.5	114 ± 6	min. 290
		<i>*non-release criteria</i>		
Product Part #		362xx.106222045		
Design and Feature		Hand specific, curved fingers, micro roughened surface and beaded cuff		
Type		Powder free, polymer coated with light chlorination and Gamma-sterilized surgical gloves		
Material		Natural Rubber Latex		
Colour		White		
Surface Treatment		Chlorination (1L-PcP)		
Powder Residue (mg/glove)		≤ 2		
Protein Content (µg/g)		≤ 50		
Single-wall Thickness (mm)		<b>Specification</b>		
Finger		0.22 ± 0.03		
Palm		0.18 ± 0.03		
*All sizes Cuff		0.15 ± 0.02		
Physical Properties		<b>Unaged</b>		<b>Aged</b>
Tensile Strength (MPa)		Min. 24		Min. 18
Ultimate Elongation (%)		Min. 750		Min. 560
Stress at 500% (MPa)		Max. 5.5		N/A
Force at break (N)		≥ 9		≥ 9
Packing Mode		1 pairs of gloves per inner wrap 1 inner wrap per pouch 50 pouch per dispenser 4 dispenser per carton		
Glove Marking		The cuff of each glove is marked with 'Hand and Size'		
Lot # Identification on Finished Goods		Lot Number Structure : YMMPPPPSS (9 Digits) Y = Year of packing (e.g. 3 for year 2023, 4 for year 2024, etc.) MM = Month of packing (e.g. 09 for Sept., 10 for Oct., etc.) PPPP = WRP's Packing Work Order# (PWO) SS = Size (55=5.5, 60=6.0, 65=6.5, 70=7.0, 75=7.5, 80=8.0, 85=8.5 and 90=9.0)		
Shelf Life		3 years upon manufacturing date		
Pre-shipment Inspection		N=13, Median; S-2, AQL 4.0 N=13, Median; S-2, AQL 4.0 G-1, AQL 0.65 N=8 N=3 G-1, AQL 1.5 G-1, AQL 2.5		
*Single-Normal Sampling Plan				
Product Conformance		<ul style="list-style-type: none"> <li>Personal Protective Equipment of Complex Design Category III, in compliance with European Regulation (EU) 2016/425, Regulation 2016/425 on personal protective equipment, as amended to apply in GB, tested to EN ISO 21420:2020, EN ISO 374-1:2016+A1:2018 Type B, EN ISO 374-2:2019, EN ISO 374-4:2019, EN ISO 374-5:2016, EN 16523-1:2015+A1:2018, CE2797</li> <li>CE Class IIa Sterile Medical Device: in compliance with 93/42/EEC</li> <li>EN455 Parts 1, 2, 3 and 4</li> <li>ASTM D3577</li> </ul>		
Quality Assurance		<ul style="list-style-type: none"> <li>MDSAP ISO 13485:2016 Quality Management System</li> <li>ISO 9001:2015 Quality Management System</li> <li>EN ISO 13485:2016</li> <li>ISO 14001:2015, Environmental Management System</li> </ul>		

WRP ASIA PACIFIC SDN. BHD.



Name : Rosnaini Binti Rashidin

Title : Product Management Manager