



E-DUO™ Tone

**Dual colored, Heavy-duty Nitrile gloves
with Raised diamond textured surface**

Key Features

- Three layers of protection**
- Raised Diamond Texture provides excellent grip**
- Black inner layer for fast detection of breaks**
- Touchscreen compatible**

Product Information

Reorder No.	E85xx-01
Material	Nitrile
Sizes	S - XXL
Length (mm)	Min. 240
Target Weight (g)	8.5
Finger Thickness Smooth Area (mm)	0.17 ± 0.03
Colour	Inner : Black Outer : Olive green
Packing Mode	50 pc/bx; 10 bx/cs <small>* Size XXL is 45 pieces/box</small>

PRODUCT CONFORMANCE

ASTM D6319
Medical Device: in compliance with MDR (EU)2017/745
(CE Class I)
EN 455: Part 1/2/3/4
CE2797

QUALITY ASSURANCE

MDSAP-ISO 13485:2016 Quality Management System
ISO 9001:2015 Quality Management System
EN ISO 13485:2016
ISO 14001:2015 Environmental Management System



WRP Asia Pacific Sdn Bhd

Lot 1, Jalan 3, Kawasan Perusahaan Bandar Baru Salak Tinggi
43900 Sepang, Selangor Darul Ehsan, Malaysia
Phone +60-3-8706-1486
Email inquiry@wrpworld.com

EC REP & Importer in the EU

Remesco Handelsges.m.b.H.
Muthgasse 36/19, A-1190 Vienna, Austria
Phone +43-1-328-50-88
Email office@remesco.com