



SPECIFICATION SHEET
E-DUO Tone

Size	Reorder No.	Dimension		
		Target weight (g/pcs)	Palm Width (mm)	Length (mm)
Small	E8501-01	7.5 ± 0.5	88 ± 5	Min. 240
Medium	E8502-01	8.5 ± 0.5	98 ± 5	Min. 240
Large	E8503-01	9.5 ± 0.5	108 ± 5	Min. 240
X-Large	E8504-01	10.0 ± 0.5	118 ± 5	Min. 240
2X-Large	E8505-01	10.5 ± 0.5	128 ± 5	Min. 240
		<i>*Non-release criteria</i>		
Design and Feature		Ambidextrous, raised diamond embossed at palm and back of hand		
Type		Powder free and non-sterile		
Material		100% Nitrile (Acrylonitrile-butadiene)		
Colour		Inner : Black Outer : Olive green (Pantone 5555 C)		
Surface Treatment		Polymer coating and chlorination		
Powder Residue (mg/glove)		≤ 2		
Single-wall Thickness (mm) <i>*All sizes</i>	Finger Palm Cuff	Smooth area		Textured area
		0.17 ± 0.03		0.23 ± 0.03
		0.15 ± 0.03		0.20 ± 0.03
		0.11 ± 0.03		0.12 ± 0.03
		<i>* Release criteria</i>		<i>* Non-release criteria</i>
Physical Properties	Tensile Strength (MPa)	Min.18		Min.14
	Elongation (%)	Min. 500		Min. 400
	Force at break (N)	≥ 6		≥ 6
Packing Mode		50 gloves by weight and 10 dispenser per carton (S-XL) 45 gloves by weight and 10 dispenser per carton (2XL)		
Glove Marking		No marking		
Lot # Identification on Finished Goods		Lot Number Structure : YMMPPPPSS (9 Digits) Y = Year of packing (e.g. 3 for year 2023, 4 for year 2024, etc.) MM = Month of packing (e.g. 09 for Sept., 10 for Oct., etc.) PPPP = WRP's Packing Work Order# (PWO) SS = Size (01 = S, 02 = M, 03 = L, 04 = XL and 05 = XXL)		
Shelf Life Claim		5 years upon manufacturing date		
Pre-shipment Inspection <i>*Single-Normal Sampling Plan</i>	Dimension	N = 13, Median; S-2, AQL 4.0		
	Physical Properties	N = 13, Median; S-2, AQL 4.0		
	1000ml Water Leak	G-1, AQL 1.5		
	Powder Residue	N = 5		
	Major Visual Inspection	G-1, AQL 2.5		
	Minor Visual Inspection	G-1, AQL 4.0		
Product Conformance		<ul style="list-style-type: none"> Medical Device: in compliance with MDR (EU) 2017/745 (CE Class II) EN ISO 21420:2020, EN ISO 374-1:2016 + A1:2018 Type B, EN ISO 374-2:2019, EN ISO 374-4:2019, EN ISO 374-5:2016, EN 16523-1:2015 + A1:2018, CE2797 EN 455 – Part 1,2,3,4 ASTM D6319 		
Quality Assurance		<ul style="list-style-type: none"> MDSAP ISO 13485:2016 Quality Management System ISO 9001:2015 Quality Management System EN ISO 13485:2016 ISO 14001:2015, Environmental Management System 		

Note: The above information is a guideline of typical performance values and characteristics of the product.

WRP ASIA PACIFIC SDN. BHD.




Name : Rosnaini Binti Rashidin
 Title : Product Management Manager