



# COMFIT®

## SURGICAL GLOVES

NR LATEX POWDER FREE

**Hand Specific**  
**Powder Free**  
**Micro Roughened Surface**

### Key Features

- Polymer coated**
- Soft and comfortable**
- Beaded cuff**

### Product Information

|                          |                      |
|--------------------------|----------------------|
| Reorder No.              | C40xx-93             |
| Material                 | Natural Rubber Latex |
| Sizes                    | 5.5 - 9.0            |
| Length (mm)              | ± 290                |
| Target Weight (Size 7.5) | 9.0g                 |
| Finger Thickness (mm)    | 0.18 ± 0.03          |
| Colour                   | Natural              |
| Packing Mode             | 50 pr/bx; 4 bx/cs    |

### PRODUCT CONFORMANCE

- ASTM D3577
- Medical Device: in compliance with 93/42/EEC
- EN 455: Part 1/2/3/4
- EN ISO 21420:2020, EN ISO 374-1:2016+A1:2018 Type B, EN ISO 374-2:2019, EN ISO 374-4:2019, EN ISO 374-5:2016, EN 16523-1:2015+A1:2018, ISO 16604:2004, CE2797

### QUALITY ASSURANCE

- MDSAP ISO 13485:2016 Quality Management System
- ISO 9001:2015 Quality Management System
- EN ISO 13485:2016
- ISO 14001:2015 Environmental Management System



### **WRP Asia Pacific Sdn Bhd**

Lot 1, Jalan 3, Kawasan Perusahaan Bandar Baru Salak Tinggi  
 43900 Sepang, Selangor Darul Ehsan, Malaysia  
 Phone +60-3-8706-1486  
 Email [inquiry@wrpworld.com](mailto:inquiry@wrpworld.com)

### **EC REP & Importer in the EU**

Remesco Handelses.m.b.H.  
 Muthgasse 36/19, A-1190 Vienna, Austria  
 Phone +43-1-328-50-88  
 Email [office@remesco.com](mailto:office@remesco.com)