



COMFIT® Premium

POWDER FREE LATEX SURGICAL GLOVES

**Hand Specific
Powder Free
Micro Roughened Surface**

Key Features

Sterilized

Beaded cuff

Polymer coated

Enhances comfort

Soft and comfortable

Enabling good dexterity

Non slick skin grip finishing

Product Information

Reorder No.	C36xx-96
Material	Natural Rubber Latex
Sizes	5.5 - 9.0
Length (mm)	± 290
Target Weight (Size 7.5)	10.5g
Finger Thickness (mm)	0.21 ± 0.03
Colour	Natural
Packing Mode	50 pr/bx; 4 bx/cs

PRODUCT CONFORMANCE

- ASTM D3577
- Medical Device: in compliance with 93/42/EEC
- EN 455: Part 1/2/3/4
- EN ISO 21420:2020, EN ISO 374-1:2016+A1:2018 Type B, EN ISO 374-2:2019, EN ISO 374-4:2019, EN ISO 374-5:2016, EN 16523-1:2015+A1:2018, ISO 16604:2004, CE2797

QUALITY ASSURANCE

- MDSAP-ISO 13485:2016 Quality Management System
- ISO 9001:2015 Quality Management System
- EN ISO 13485:2016
- ISO 14001:2015 Environmental Management System



WRP Asia Pacific Sdn Bhd

Lot 1, Jalan 3, Kawasan Perusahaan Bandar Baru Salak Tinggi
43900 Sepang, Selangor Darul Ehsan, Malaysia
Phone +60-3-8706-1486
Email inquiry@wrpworld.com

EC REP & Importer in the EU

Remesco Handelsges.m.b.H.
Muthgasse 36/19, A-1190 Vienna, Austria
Phone +43-1-328-50-88
Email office@remesco.com