



EPIC^{MT}

SURGICAL GLOVES

NR LATEX • POWDER FREE • MICRO TEXTURED GRIP

Hand Specific
Powder Free
Micro Roughened Surface

Key Features

- Sterilized**
- Beaded cuff**
- Polymer coated**
- Enformance performance**
- Superior tactile sensitivity**

Product Information

Reorder No.	EP-12 -xx
Material	Natural Rubber Latex
Sizes	5.5 - 9.0
Length (mm)	± 290
Target Weight (Size 7.5)	11.5g
Finger Thickness (mm)	0.22 ± 0.03
Colour	White
Packing Mode	50 pr/bx; 6 bx/cs 50 pr/bx; 4 bx/cs

PRODUCT CONFORMANCE

- ASTM D3577
- Medical Device: in compliance with 93/42/EEC
- EN 455: Part 1/2/3/4
- EN ISO 21420:2020, EN ISO 374-1:2016+A1:2018 Type B, EN ISO 374-2:2019, EN ISO 374-4:2019, EN ISO 374-5:2016, EN 16523-1:2015+A1:2018, ISO 16604:2004, CE2797

QUALITY ASSURANCE

- MDSAP-ISO 13485:2016 Quality Management System
- ISO 9001:2015 Quality Management System
- EN ISO 13485:2016
- ISO 14001:2015 Environmental Management System



WRP Asia Pacific Sdn Bhd

Lot 1, Jalan 3, Kawasan Perusahaan Bandar Baru Salak Tinggi
 43900 Sepang, Selangor Darul Ehsan, Malaysia
 Phone +60-3-8706-1486
 Email inquiry@wrpworld.com

EC REP & Importer in the EU

Remesco Handelsges.m.b.H.
 Muthgasse 36/19, A-1190 Vienna, Austria
 Phone +43-1-328-50-88
 Email office@remesco.com