



PROFEEL®

DHD™ Platinum
Powder Free Latex Surgical Gloves

Hand Specific
Powder Free
Damp-Hand-Donnable (DHD™)
Micro Roughened Surface

Key Features

- Sterilized**
- Beaded cuff**
- Polymer coated**
- Non-slick skin grip**
- Soft and comfortable**
- Superior tactile sensitivity**
- Chemotherapy tested**

Product Information

Reorder No.	P33xx-01
Material	Natural Rubber Latex
Sizes	5.5 - 9.0
Length (mm)	± 290
Target Weight (Size 7.5)	10.8g
Finger Thickness (mm)	0.21 ± 0.03
Colour	Natural
Packing Mode	50 pr/bx; 4 bx/cs

PRODUCT CONFORMANCE

- ASTM D3577
- Meets ASTM D6978 for Chemotherapy testing
- Medical Device: in compliance with 93/42/EEC
- EN 455: Part 1/2/3/4
- EN ISO 21420:2020, EN ISO 374-1:2016+A1:2018 Type B, EN ISO 374-2:2019, EN ISO 374-4:2019, EN ISO 374-5:2016, EN 16523-1:2015+A1:2018, ISO 16604:2004, CE2797
- Health Canada Compliance

QUALITY ASSURANCE

- MDSAP-ISO 13485:2016 Quality Management System
- ISO 9001:2015 Quality Management System
- EN ISO 13485:2016
- ISO 14001:2015 Environmental Management System



WRP Asia Pacific Sdn Bhd 1 4 7 8 1 7 V
Lot 1, Jalan 3, Kawasan Perusahaan Bandar Baru Salak Tinggi,
43900 Sepang, Selangor Darul Ehsan, MALAYSIA
Office +60-3-8706-1486
Email inquiry@wrpworld.com

EC REP & Importer in the EU
REMESCO Handelsges.m.b.H.
Muthgasse 36/19, A-1190 Vienna, AUSTRIA
Office +43-1-328-50-88
Fax +43-1-328-29-49
Email office@remesco.com