



# PROFEEL®

**DHD™ Polyisoprene Sensitive**  
Powder Free Non-Latex Surgical Gloves

**Hand Specific  
Powder Free  
Damp-Hand-Donnable (DHD™)**

## Key Features

- Sterilized**
- Latex free**
- Beaded cuff**
- Smooth finish**
- Polymer coated**
- Optimum sensitivity**
- Optimal instrument handling and grip**

## Product Information

Reorder No.	P33xx-57
Material	Polyisoprene
Sizes	5.5 - 9.0
Length (mm)	± 290
Target Weight (Size 7.5)	11.8g
Finger Thickness (mm)	0.20 ± 0.02
Colour	Green
Packing Mode	50 pr/bx; 4 bx/cs

## PRODUCT CONFORMANCE

- ASTM D3577
- Medical Device: in compliance with 93/42/EEC
- EN 455: Part 1/2/3/4
- EN ISO 21420:2020, EN ISO 374-1:2016+A1:2018 Type B, EN ISO 374-2:2019, EN ISO 374-4:2019, EN ISO 374-5:2016, EN 16523-1:2015+A1:2018, ISO 16604:2004, CE2797
- Health Canada Compliance

## QUALITY ASSURANCE

- MDSAP ISO 13485:2016 Quality Management System
- ISO 9001:2015 Quality Management System
- EN ISO 13485:2016
- ISO 14001:2015 Environmental Management System



**WRP Asia Pacific Sdn Bhd** 1 4 7 8 1 7 V  
Lot 1, Jalan 3, Kawasan Perusahaan Bandar Baru Salak Tinggi,  
43900 Sepang, Selangor Darul Ehsan, MALAYSIA  
Office +60-3-8706-1486  
Email [inquiry@wrpworld.com](mailto:inquiry@wrpworld.com)

**EC REP** & Importer in the EU  
**REMESCO Handelsges.m.b.H.**  
Muthgasse 36/19, A-1190 Vienna, AUSTRIA  
Office +43-1-328-50-88  
Fax +43-1-328-29-49  
Email [office@remesco.com](mailto:office@remesco.com)