



PROFEEL®

DOUBLE GLOVING SYSTEM

FOR PUNCTURE DETECTION
DHD™ LATEX SURGICAL UNDERGLOVE AND OVERGLOVE POWDER FREE

**Hand Specific
Powder Free
Damp-Hand-Donnable (DHD™)**

Key Features

Sterilized

Beaded cuff

Polymer coated

Micro roughened surface

Standard thickness / Micro thickness

Product Information

	<u>Overglove</u>	<u>Underglove</u>
Reorder No.	P33xx-61	P33xx-61
Material	Natural Rubber Latex	Natural Rubber Latex
Sizes	5.5 - 9.0	5.5 - 9.0
Length (mm)	± 290	± 290
Target Weight (Size 7.5)	10.8g	9.9g
Finger Thickness (mm)	0.21 ± 0.03	0.20 ± 0.03
Colour	Natural	Dark Green
Packing Mode	25pr/bx; 4 bx/cs	25pr/bx; 4 bx/cs

PRODUCT CONFORMANCE

- ASTM D3577
- Medical Device: in compliance with 93/42/EEC
- EN 455: Part 1/2/3/4
- EN ISO 21420:2020, EN ISO 374-1:2016+A1:2018 Type B, EN ISO 374-2:2019, EN ISO 374-4:2019, EN ISO 374-5:2016, EN 16523-1:2015+A1:2018, ISO 16604:2004, CE2797

QUALITY ASSURANCE

- MDSAP-ISO 13485:2016 Quality Management System
- ISO 9001:2015 Quality Management System
- EN ISO 13485:2016
- ISO 14001:2015 Environmental Management System



WRP Asia Pacific Sdn Bhd

Lot 1, Jalan 3, Kawasan Perusahaan Bandar Baru Salak Tinggi
43900 Sepang, Selangor Darul Ehsan, Malaysia
Phone +60-3-8706-1486
Email inquiry@wrpworld.com

EC REP & Importer in the EU

Remesco Handelsges.m.b.H.
Muthgasse 36/19, A-1190 Vienna, Austria
Phone +43-1-328-50-88
Email office@remesco.com