



RadiaXon®

Radiation Attenuating Surgical Gloves

LEAD FREE NR Latex Powder Free Gloves

NR Latex
Hand Specific
Powder Free

Key Features

Sterilized

Disposable

Beaded cuff

Lead free and non toxic

Environmentally friendly

Product Information

Reorder No.	R32xx-28
Material	Natural Rubber Latex
Sizes	6.0 - 9.0
Length (mm)	± 290
Target Weight (Size 7.5)	47g
Finger Thickness (mm)	0.33 ± 0.03
Colour	Brown
Packing Mode	5 pr/bx; 4 bx/cs

PRODUCT CONFORMANCE

- ASTM D3577
- Medical Device: in compliance with 93/42/EEC
- EN 455: Part 1/2/3/4
- EN ISO 21420:2020, EN 421:2010, EN ISO 374-1:2016+A1:2018 Type B, EN ISO 374-2:2019, EN ISO 374-4:2019, EN ISO 374-5:2016, EN 16523-1:2015+A1:2018, ISO 16604:2004, EN 61331-3:2014, CE2797

QUALITY ASSURANCE

- MDSAP-ISO 13485:2016 Quality Management System
- ISO 9001:2015 Quality Management System
- EN ISO 13485:2016
- ISO 14001:2015 Environmental Management System



WRP Asia Pacific Sdn Bhd

Lot 1, Jalan 3, Kawasan Perusahaan Bandar Baru Salak Tinggi
43900 Sepang, Selangor Darul Ehsan, Malaysia
Phone +60-3-8706-1486
Email inquiry@wrpworld.com

EC REP & Importer in the EU

Remesco Handelsges.m.b.H.
Muthgasse 36/19, A-1190 Vienna, Austria
Phone +43-1-328-50-88
Email office@remesco.com