



RadiaXon®

PI Radiation Attenuating Surgical Gloves

LEAD FREE • NON-LATEX Polyisoprene Powder Free Gloves

**Polyisoprene
Hand Specific
Powder Free**

Key Features

Sterilized

Beadless cuff

Latex free

Disposable

Soft and comfortable

Lead free and non toxic

Environmentally friendly

Product Information

Reorder No.	R35xx-85
Material	Polyisoprene
Sizes	6.0 - 9.0
Length (mm)	± 290
Target Weight (Size 7.5)	47g
Finger Thickness (mm)	0.33 ± 0.03
Colour	Black
Packing Mode	5 pr/bx; 4 bx/cs

PRODUCT CONFORMANCE

- ASTM D3577
- Medical Device: in compliance with 93/42/EEC
- EN 455: Part 1/2/3/4
- EN ISO 21420:2020, EN 421:2010, EN ISO 374-1:2016+A1:2018 Type B, EN ISO 374-2:2019, EN ISO 374-4:2019, EN ISO 374-5:2016, EN 16523-1:2015+A1:2018, ISO 16604:2004, EN 61331-3:2014, CE2797

QUALITY ASSURANCE

- MDSAP-ISO 13485:2016 Quality Management System
- ISO 9001:2015 Quality Management System
- EN ISO 13485:2016
- ISO 14001:2015 Environmental Management System



WRP Asia Pacific Sdn Bhd

Lot 1, Jalan 3, Kawasan Perusahaan Bandar Baru Salak Tinggi
43900 Sepang, Selangor Darul Ehsan, Malaysia
Phone +60-3-8706-1486
Email inquiry@wrpworld.com

EC REP & Importer in the EU

Remesco Handelsges.m.b.H.
Muthgasse 36/19, A-1190 Vienna, Austria
Phone +43-1-328-50-88
Email office@remesco.com