



PROFEEL®

DHD™ Platinum Plus

Powder Free Latex Surgical Gloves

Damp-Hand-Donnable (DHD™)
Micro-roughened surface
Powder Free

Key Features

- Chemotherapy and fentanyl tested**
- Polymer coated**
- Skin-grip finish**
- Soft and comfortable**
- Superior tactile sensitivity**
- Beaded cuff**

Product Information

Reorder No.	P40xx-02
Material	Natural Rubber Latex
Sizes	5.5 - 9.0
Length (mm)	290
Finger Thickness (mm)	0.21
Colour	Natural
Packing Mode	50 pr/bx; 4 bx/cs

PRODUCT CONFORMANCE

- ASTM D3577
- Meets ASTM D6978 for Chemotherapy and Fentanyl testing
- Medical Device: in compliance with MDR (EU) 2017/745
- UK Medical Device Regulations 2002 Part II (UKCA Class IIa)
- EN 455: Part 1/2/3/4
- EN ISO 21420:2020, EN ISO 374-1:2016+A1:2018 Type B, EN ISO 374-2:2019, EN ISO 374-4:2019, EN ISO 374-5:2016, EN 16523-1:2015+A1:2018, ISO 16604:2004, CE2797.

QUALITY ASSURANCE

- MDSAP ISO 13485:2016 Quality Management System
- ISO 9001:2015 Quality Management System
- EN ISO 13485:2016
- ISO 14001:2015 Environmental Management System



WRP Asia Pacific Sdn Bhd
 Lot 1, Jalan 3, Kawasan Perusahaan Bandar Baru Salak Tinggi
 43900 Sepang, Selangor Darul Ehsan, Malaysia
 Phone: +60 3 8706 1486
 Email: inquiry@wrpworld.com

[EC] REP & Importer in the EU
 Remesco Handelsges.m.b.H.
 Muthgasse 36/19, A-1190 Vienna, Austria
 Phone: +43 1328 5088
 Email: office@remesco.com

[UK] REP
 Obelis UK Ltd
 Sandford Gate, Oxford OX4 6LB, UK
 Phone: +32 (0) 2732 5954
 Email: mail@obelis.net

[CH] REP
 Obelis Swiss G.m.b.H.
 Ruessenstrasse 12, 6340 Baar / ZG, CH
 Phone: +41 41544 1526
 Email: info@obelis.ch