



**NR Latex  
Powder Free**

### Key Features

- Lead free and non toxic**
- Environmentally friendly**
- Disposable**
- Beaded cuff**

### Product Information

|                       |                      |
|-----------------------|----------------------|
| Reorder No.           | XR49xx-00            |
| Material              | Natural Rubber Latex |
| Sizes                 | 6.0 - 9.0            |
| Length (mm)           | 290                  |
| Finger Thickness (mm) | 0.33                 |
| Colour                | Khaki                |
| Packing Mode          | 5 pr/bx; 4 bx/cs     |

### **PRODUCT CONFORMANCE**

- ASTM D7866-14a (Primary Radiation)
- ASTM F3094-14 (Scattered Radiation)

### **QUALITY ASSURANCE**

- MDSAP-ISO 13485:2016 Quality Management System
- ISO 9001:2015 Quality Management System
- EN ISO 13485:2016
- ISO 14001:2015 Environmental Management System



SINGLE USE  
ONLY



STERILISED BY  
IRRADIATION

 **WRP Asia Pacific Sdn Bhd**

Lot 1, Jalan 3, Kawasan Perusahaan Bandar Baru Salak Tinggi  
43900 Sepang, Selangor Darul Ehsan, Malaysia  
Phone: +60 3 8706 1486  
Email: [inquiry@wrpworld.com](mailto:inquiry@wrpworld.com)