



**COMFIT®**  
SURGICAL GLOVES  
NR LATEX POWDER FREE

**Damp-Hand-Donnable (DHD™)**  
**Powder Free**  
**Micro Roughened Surface**

**Key Features**

- Polymer coated**
- Soft and comfortable**
- Beaded cuff**

**Product Information**


|                       |                      |
|-----------------------|----------------------|
| Reorder No.           | C40xx-93             |
| Material              | Natural Rubber Latex |
| Sizes                 | 5.5 - 9.0            |
| Length (mm)           | 290                  |
| Finger Thickness (mm) | 0.18                 |
| Colour                | Natural              |
| Packing Mode          | 50 prs/bx; 4 bx/cs   |

**PRODUCT CONFORMANCE**

- ASTM D3577
- ASTM D6978 - Tested for use with Chemotherapy Drugs & Fentanyl Citrate
- Medical Device: in compliance with MDR (EU) 2017/745 (CE Class IIa)
- UK Medical Device Regulations 2002 Part II (UKCA Class IIa)
- EN 455: Part 1/2/3/4
- EN ISO 21420:2020, EN ISO 374-1:2016+A1:2018 Type B, EN ISO 374-2:2019, EN ISO 374-4:2019, EN ISO 374-5:2016, EN 16523-1:2015+A1:2018, ISO 16604:2004 & CE2797
- Health Canada Compliance
- EUDR Compliant Latex available upon request
- FDA 510(k) available

**QUALITY ASSURANCE**

- MDSAP ISO 13485:2016 Quality Management System
- ISO 9001:2015 Quality Management System
- EN ISO 13485:2016 Quality Management System - Regulatory Purpose
- ISO 14001:2015 Environmental Management System



**WWRP Asia Pacific Sdn Bhd**  
Lot 1, Jalan 3, Kawasan Perusahaan Bandar Baru Salak Tinggi  
43900 Sepang, Selangor Darul Ehsan, Malaysia  
Phone: +603 8706 1486  
Email: inquiry@wrpworld.com

**EC REP** & Importer in the EU  
**Remesco Handelsges.m.b.H.**  
Muthgasse 36/19, A-1190 Vienna, Austria  
Phone: +43 1328 5088  
Email: office@remesco.com

**UK REP**  
**Obelis UK Ltd**  
Sandford Gate, Oxford OX4 6LB, UK  
Phone: +32 (0) 2732 5954  
Email: mail@obelis.net

**CH REP**  
**Obelis Swiss G.m.b.H.**  
Ruessenstrasse 12, 6340 Baar / ZG, CH  
Phone: +41 41544 1526  
Email: info@obelis.ch