



COMFIT® SURGICAL GLOVES NR LATEX POWDERED

***Powdered
Micro Roughened Surface***

Key Features

Superior sensitivity

Beaded cuff

Product Information

Reorder No.	C33xx-10
Material	Natural Rubber Latex
Sizes	5.5 - 9.0
Length (mm)	290
Finger Thickness (mm)	0.18
Colour	Natural
Packing Mode	50 prs/bx; 4 bx/cs

PRODUCT CONFORMANCE

- EN 455: Part 1/2/3/4
- EN ISO 21420:2020, EN ISO 374-1:2016+A1:2018 Type B, EN ISO 374-2:2019, EN ISO 374-4:2019, EN ISO 374-5:2016, EN 16523-1:2015+A1:2018, ISO 16604:2004, CE2797.
- EUDR Compliant Latex available upon request.
- Powdered Surgical Gloves are available for regions where CE MDR and Health Canada certification are not required. The product remains certified under CE PPE Category III, UKCA and MDA.

QUALITY ASSURANCE

- MDSAP ISO 13485:2016 Quality Management System
- ISO 9001:2015 Quality Management System
- EN ISO 13485:2016 Quality Management System - Regulatory Purpose
- ISO 14001:2015 Environmental Management System



WRP Asia Pacific Sdn Bhd

Lot 1, Jalan 3, Kawasan Perusahaan Bandar Baru Salak Tinggi
43900 Sepang, Selangor Darul Ehsan, Malaysia
Phone: +603 8706 1486
Email: inquiry@wrpworld.com