



DERMAGRIP®

X-Black Nitrile Examination Gloves

Powder Free
Nitrile
Ambidextrous
Fully Textured Surface

Key Features

Chemotherapy & Fentanyl tested

Soft comfortable formulation

Excellent dexterity and elasticity

Extra strong and durable

Beaded cuff

Product Information


Reorder No.	D15xx-93
Material	Nitrile
Sizes	S - XXL
Length (mm)	240
Finger Thickness (mm)	0.18
Colour	Black
Packing Mode	100 pcs/bx, 10 bx/cs

PRODUCT CONFORMANCE

- ASTM D6319
- ASTM D6978 - Tested for use with Chemotherapy Drugs & Fentanyl Citrate
- Medical Device: in compliance with MDR (EU)2017/745 (CE Class I)
- UK Medical Device Regulations 2002 Part II (UKCA Class I)
- EN 455: Part 1/2/3/4
- EN ISO 21420:2020, EN ISO 374-1:2016+A1:2018 Type B, EN ISO 374-2:2019, EN ISO 374-4:2019, EN ISO 374-5:2016, EN 16523-1:2015+A1:2018, ISO 16604:2004 & CE 2797
- Food Contact Compliance FDA 21 CFR 177.2600 & EU ResAP (2004) 4
- FDA 510(k) available

QUALITY ASSURANCE

- MDSAP ISO 13485:2016 Quality Management System
- ISO 9001:2015 Quality Management System
- EN ISO 13485:2016 Quality Management System - Regulatory Purpose
- ISO 14001:2015 Environmental Management System



WRP Asia Pacific Sdn Bhd
Lot 1, Jalan 3, Kawasan Perusahaan Bandar Baru Salak Tinggi
43900 Sepang, Selangor Darul Ehsan, Malaysia
Phone: +603 8706 1486
Email: inquiry@wrpworld.com

[EC REP] & Importer in the EU
Remesco Handelsges.m.b.H.
Muthgasse 36/19, A-1190 Vienna, Austria
Phone: +43 1328 5088
Email: office@remesco.com

[UK REP]
Obelis UK Ltd
Sandford Gate, Oxford OX4 6LB, UK
Phone: +32 (0) 2732 5954
Email: mail@obelis.net

[CH REP]
Obelis Swiss G.m.b.H.
Ruessenstrasse 12, 6340 Baar / ZG, CH
Phone: +41 41544 1526
Email: info@obelis.ch