



DERMAGRIP®
HIGH RISK EXAMINATION GLOVES
NR LATEX POWDER FREE

NR Latex
Ambidextrous
Powder Free
Finger Textured Surface

Key Features

- No unpleasant odor
- Chemotherapy tested
- Extra-long and non-sterile protection
- Thicker and stronger
- Protection in high risk situations
- Beaded cuff

Product Information


Reorder No.	D14xx-14
Material	Natural Rubber Latex
Sizes	S - XL
Length (mm)	290
Finger Thickness (mm)	0.36
Colour	Blue
Packing Mode	50 pcs/bx, 10 bx/cs

PRODUCT CONFORMANCE

- ASTM D3578
- ASTM D6978 - Tested for use with Chemotherapy Drugs
- Medical Device: in compliance with MDR (EU) 2017/745 (CE Class I)
- UK Medical Device Regulations 2002 Part II (UKCA Class I)
- EN 455: Part 1/2/3/4
- EN ISO 21420:2020, EN ISO 374-1:2016+A1:2018 Type A, EN ISO 374-2:2019, EN ISO 374-4:2019, EN ISO 374-5:2016, EN 16523-1:2015+A1:2018, ISO 16604:2004 & CE2797
- Food Contact Compliance: EU No. 10/2011
- Health Canada Compliance
- EUDR Compliant Latex available upon request
- FDA 510(k) available

QUALITY ASSURANCE

- MDSAP ISO 13485:2016 Quality Management System
- ISO 9001:2015 Quality Management System
- EN ISO 13485:2016 Quality Management System - Regulatory Purpose
- ISO 14001:2015 Environmental Management System



WRP Asia Pacific Sdn Bhd
Lot 1, Jalan 3, Kawasan Perusahaan Bandar Baru Salak Tinggi
43900 Sepang, Selangor Darul Ehsan, Malaysia
Phone: +603 8706 1486
Email: inquiry@wrpworld.com

EC REP & Importer in the EU
Remesco Handelsges.m.b.H.
Muthgasse 36/19, A-1190 Vienna, Austria
Phone: +43 1328 5088
Email: office@remesco.com

UK REP
Obelis UK Ltd
Sandford Gate, Oxford OX4 6LB, UK
Phone: +32 (0) 2732 5954
Email: mail@obelis.net

CH REP
Obelis Swiss G.m.b.H.
Ruessenstrasse 12, 6340 Baar / ZG, CH
Phone: +41 41544 1526
Email: info@obelis.ch