



DERMAGRIP® DuoTone

**Dual colored, Extra-strength
Nitrile gloves with Textured fingers**

Key Features

Chemotherapy & Fentanyl tested

Three layers of protection

25% longer than standard exam gloves

Hi-Viz inner layer for fast detection of breaks

Textured fingers provide excellent grip

Touchscreen compatible

Product Information

Reorder No.	D14xx-03
Material	Nitrile
Sizes	S - XXL
Length (mm)	280
Finger Thickness (mm)	0.20
Colour	Inner : Hi-Viz Green Outer : White
Packing Mode	50 pcs/bx; 10 bx/cs * Size XL and XXL is 45 pieces/box

PRODUCT CONFORMANCE

- ASTM D6319
- ASTM D6978 - Tested for use with Chemotherapy Drugs & Fentanyl Citrate
- Medical Device: in compliance with MDR (EU)2017/745 (CE Class I)
- UK Medical Device Regulations 2002 Part II (UKCA Class I)
- EN 455: Part 1/2/3/4
- EN ISO 21420:2020, EN ISO 374-1:2016+A1:2018 Type B, EN ISO 374-2:2019, EN ISO 374-4:2019, EN ISO 374-5:2016, EN 16523-1:2015+A1:2018, ISO 16604:2004 & CE 2797
- Food Contact Compliance FDA 21 CFR 177.2600
- Health Canada Compliance
- FDA 510(k) available

QUALITY ASSURANCE

- MDSAP ISO 13485:2016 Quality Management System
- ISO 9001:2015 Quality Management System
- EN ISO 13485:2016 Quality Management System - Regulatory Purpose
- ISO 14001:2015 Environmental Management System



WRP Asia Pacific Sdn Bhd
Lot 1, Jalan 3, Kawasan Perusahaan Bandar Baru Salak Tinggi
43900 Sepang, Selangor Darul Ehsan, Malaysia
Phone: +603 8706 1486
Email: inquiry@wrpworld.com

EC REP & Importer in the EU
Remesco Handelsges.m.b.H.
Muthgasse 36/19, A-1190 Vienna, Austria
Phone: +43 1328 5088
Email: office@remesco.com

UK REP
Obelis UK Ltd
Sandford Gate, Oxford OX4 6LB, UK
Phone: +32 (0) 2732 5954
Email: mail@obelis.net

CH REP
Obelis Swiss G.m.b.H.
Ruessenstrasse 12, 6340 Baar / ZG, CH
Phone: +41 41544 1526
Email: info@obelis.ch