



RadiaXon®

Radiation Attenuating Surgical Gloves

LEAD FREE NR Latex Powder Free Gloves

NR Latex
Powder Free

Key Features

Disposable

Lead free and non toxic

Environmentally friendly

Beaded cuff

Product Information


Reorder No.	R32xx-28
Material	Natural Rubber Latex
Sizes	6.0 - 9.0
Length (mm)	290
Finger Thickness (mm)	0.33
Colour	Khaki
Packing Mode	5 prs/bx; 4 bx/cs

PRODUCT CONFORMANCE

- Medical Device: in compliance with MDR (EU) 2017/745 (CE Class IIa)
- UK Medical Device Regulations 2002 Part II (UKCA Class IIa)
- EN 455: Part 1/2/3/4
- EN ISO 21420:2020, EN ISO 374-1:2016+A1:2018 Type B, EN ISO 374-2:2019, EN ISO 374-4:2019, EN ISO 374-5:2016, EN 16523-1:2015+A1:2018, ISO 16604:2004, EN 61331-3:2014, EN 421:2010 & CE2797
- EUDR Compliant Latex available upon request.

QUALITY ASSURANCE

- MDSAP ISO 13485:2016 Quality Management System
- ISO 9001:2015 Quality Management System
- EN ISO 13485:2016 Quality Management System - Regulatory Purpose
- ISO 14001:2015 Environmental Management System



WRP Asia Pacific Sdn Bhd
Lot 1, Jalan 3, Kawasan Perusahaan Bandar Baru Salak Tinggi
43900 Sepang, Selangor Darul Ehsan, Malaysia
Phone: +603 8706 1486
Email: inquiry@wrpworld.com

EC REP & Importer in the EU
Remesco Handelsges.m.b.H.
Muthgasse 36/19, A-1190 Vienna, Austria
Phone: +43 1328 5088
Email: office@remesco.com

UK REP
Obelis UK Ltd
Sandford Gate, Oxford OX4 6LB, UK
Phone: +32 (0) 2732 5954
Email: mail@obelis.net

CH REP
Obelis Swiss G.m.b.H.
Ruessenstrasse 12, 6340 Baar / ZG, CH
Phone: +41 41544 1526
Email: info@obelis.ch