



RadiaXon®

PI Radiation Attenuating Surgical Gloves

LEAD FREE • NON-LATEX Polyisoprene Powder Free Gloves

Polyisoprene
Powder Free

Key Features

- Latex free
- Disposable
- Soft and comfortable
- Lead free and non toxic
- Enviromentally friendly
- Beadless cuff

Product Information

Reorder No.	R35xx-85
Material	Polyisoprene
Sizes	6.0 - 9.0
Length (mm)	290
Finger Thickness (mm)	0.33
Colour	Black
Packing Mode	5 prs/bx; 4 bx/cs

PRODUCT CONFORMANCE

- Medical Device: in compliance with MDR (EU) 2017/745 (CE Class IIa)
- UK Medical Device Regulations 2002 Part II (UKCA Class IIa)
- EN 455: Part 1/2/3/4
- EN ISO 21420:2020, EN ISO 374-1:2016+A1:2018 Type B, EN ISO 374-2:2019, EN ISO 374-4:2019, EN ISO 374-5:2016, EN 16523-1:2015+A1:2018, ISO 16604:2004, EN 61331-3:2014, EN 421:2010 & CE2797

QUALITY ASSURANCE

- MDSAP ISO 13485:2016 Quality Management System
- ISO 9001:2015 Quality Management System
- EN ISO 13485:2016 Quality Management System - Regulatory Purpose
- ISO 14001:2015 Environmental Management System


MEDICAL DEVICE


Medical Device
REGULATION
(EU) 2017/745


Personal Protective
Equipment
REGULATION
(EU) 2016/425


2797


0086


ISO 374-1/1 TYPE B


ISO 374-5:2016


VIRUS


USER INFORMATION



 **WRP Asia Pacific Sdn Bhd**
Lot 1, Jalan 3, Kawasan Perusahaan Bandar Baru Salak Tinggi
43900 Sepang, Selangor Darul Ehsan, Malaysia
Phone: +603 8706 1486
Email: inquiry@wrpworld.com

 & Importer in the EU
Remesco Handelsges.m.b.H.
Muthgasse 36/19, A-1190 Vienna, Austria
Phone: +43 1328 5088
Email: office@remesco.com


Obelis UK Ltd
Sandford Gate, Oxford OX4 6LB, UK
Phone: +32 (0) 2732 5954
Email: mail@obelis.net


Obelis Swiss G.m.b.H.
Ruessenstrasse 12, 6340 Baar / ZG, CH
Phone: +41 41544 1526
Email: info@obelis.ch